



	Technical Data Sheet	
Product	Anti-HSV-1/HSV-2 glycoprotein B Purified	
Cat. Number/Size	11-330-C025	0.025 mg
	11-330-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	HSV-1/HSV-2 glycoprotein B	
Clone	T111	
Format	Purified	
Reactivity	Herpes Simplex Virus	
Application	ICC, ELISA	
Isotype	Mouse IgG2b	
Specificity	The antibody T111 recognizes gB antigen of both HSV1 and HSV2 (Herpes Simplex Virus type 1 and 2), a dsDNA virus, member of Herpesviridae family.	
Other names	HSVGB, HSV-1 gB, HSV-2 gB	
Immunogen	Extract of HSV-1 infected VERO (green monkey kidney) cells.	
Entrez Gene ID	24271469	
Gene name	UL27	
NCBI Full Gene Name	envelope glycoprotein B	
UniProt ID	P10211	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.