



	Technical Data Sheet	
Product	Anti-DAXX Purified	
Cat. Number/Size	11-517-C025	0.025 mg
	11-517-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	DAXX	
Clone	DAXX-03	
Format	Purified	
Reactivity	Human	
Application	IP, WB, ICC	
Application details	Western blotting: Recommended dilution: 2 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody DAXX-03 reacts with C-terminal part of Daxx, a 110 kDa intracellular protein.	
Other names	DAP6, EAP1, BING2	
Immunogen	His6-tagged protein corresponding to amino acids 580-740 of human Daxx.	
Entrez Gene ID	1616	
Gene name	DAXX	
NCBI Full Gene Name	death domain associated protein	
UniProt ID	Q9UER7	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.