



	Technical Data Sheet	
Product	Anti-Hu CD47 Purified	
Cat. Number/Size	11-225-C025	0.025 mg
	11-225-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD47	
Clone	MEM-122	
Format	Purified	
Reactivity	Non-human primates, Human, Pig	
Application	FC (QC tested), WB, IHC-F	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml. Western blotting: Non-reducing conditions.	
Isotype	Mouse IgM	
Specificity	The antibody MEM-122 reacts with an extracellular epitope of CD47 (Integrin Associated Protein), a 50-55 kDa membrane adhesion molecule (thrombospondin receptor; immunoglobulin supergene family) expressed on leukocytes, platelets and erythrocytes. It is also expressed on epithelial cells, endothelial cells, fibroblasts and many tumor cell lines.	
Other names	IAP, OA3, MER6	
Workshop	HLDA VI: WS Code AS A051	
Immunogen	COS-7 (African green monkey) cells	
Entrez Gene ID	961	
Gene name	CD47	
NCBI Full Gene Name	CD47 molecule	
UniProt ID	Q08722	
Concentration	1 mg/ml	
Preparation	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).	
Formulation	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.