



	Technical Data Sheet	
Product	Anti-N-sulfated heparan sulfate Biotin	
Cat. Number/Size	1B-821-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	N-sulfated heparan sulfate	
Clone	HepSS-1	
Format	Biotin	
Reactivity	Species independent	
Application	WB, ICC, FC (QC tested), IHC, IP	
Application details	Flow cytometry: Recommended dilution: 3-12 µg/ml.	
Isotype	Mouse IgM	
Specificity	The mouse monoclonal antibody HepSS-1 (also known as HepSS1) recognizes N-sulfated heparan sulfate (extracellular antigen) present in many species.	
Other names	N-sulfo-rich heparan sulfate	
Immunogen	MethA murine fibrosarcoma	
Concentration	1 mg/ml	
Preparation	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.	
Formulation	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.