

Technical Data Sheet

Product	<b>Anti-Hu CD178 Purified Low Endotoxin</b>	
Cat. Number/Size	<b>12-754-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	CD178
Clone	NOK-1
Format	Purified Low Endotoxin
Reactivity	Human
Application	FC (QC tested), IP, WB, ICC, FUNC
Application details	Functional application: Blocking. Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody NOK-1 recognizes an extracellular epitope of CD178 / Fas-L, an approximately 40 kDa transmembrane glycoprotein expressed on neutrophils, monocytes, and activated T and NK cells.
Other names	APTL, FASL, FAS-ligand, CD95L, TNFSF6, TNLG1A
Workshop	HLDA VII: 70322
Immunogen	L5178Y mouse T lymphoma cells expressing recombinant human CD178
Entrez Gene ID	356
Gene name	FASLG
NCBI Full Gene Name	Fas ligand
UniProt ID	P48023

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.