

Technical Data Sheet

Product

Anti-Notch1 Purified

Cat. Number/Size

11-651-C025

0.025 mg

11-651-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	Notch1
Clone	mN1A
Format	Purified
Reactivity	Mouse, Human
Negative species	Rat
Application	FC (QC tested), IP, WB, IHC(F), ICC
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody mN1A recognizes intracellular domain of Notch 1 protein, mainly its activated form. The unprocessed Notch 1 protein is recognized with lower affinity.
Other names	AOS5, TAN1, hN1, AOVD1
Immunogen	GST fusion protein containing cdc10-NCR region of mouse Notch1
Entrez Gene ID	4851
Gene name	NOTCH1
NCBI Full Gene Name	notch receptor 1
UniProt ID	P46531

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.