



	Technical Data Sheet	
Product	Anti-Hu CD10 Purified	
Cat. Number/Size	11-413-C025	0.025 mg
	11-413-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD10	
Clone	LT10	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The antibody LT10 reacts with an extracellular epitope of CD10 (CALLA - Common acute lymphatic leukemia antigen), a 100 kDa type II integral membrane protein.	
Other names	CALLA, Neprilysin, Neutral endopeptidase, Enkephalinase, Atriopeptidase, MME	
Immunogen	mouse NALM-6 leukemia pre-B cell line (tissue/cell preparation)	
Entrez Gene ID	4311	
Gene name	MME	
NCBI Full Gene Name	membrane metalloendopeptidase	
UniProt ID	P08473	
Concentration	1 mg/ml	
Preparation	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.