

Technical Data Sheet

Product

**Anti-Hu Ig Lambda Light Chain Purified**

Cat. Number/Size

**11-505-C025**

**0.025 mg**

**11-505-C100**

**0.1 mg**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

Antigen

Ig Lambda Light Chain

Clone

4C2

Format

Purified

Reactivity

Human

Negative species

Rabbit, Hamster, Sheep, Guinea pig, Goat

Application

FC (QC tested), IHC(P), IHC(F), ELISA

Application details

Flow cytometry: Recommended dilution: 2-8 µg/ml; positive control: RAMOS human cell line.

Isotype

Mouse IgG1

Specificity

The antibody 4C2 reacts with lambda light chains (22.5 kDa) of human immunoglobulin.

Other names

Immunoglobulin lambda, Igl

Concentration

1 mg/ml

Preparation

Purified by protein-A affinity chromatography.

Formulation

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage and handling

Store at 2-8°C. Do not freeze.

Images and References

[www.exbio.cz](http://www.exbio.cz)

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.