

Technical Data Sheet

Product	Anti-Hu CD26 Purified	
Cat. Number/Size	11-608-C025	0.025 mg
	11-608-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD26
Clone	BA5b
Format	Purified
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-2 µg/ml; positive control: human peripheral blood leukocytes.
Isotype	Mouse IgG2a kappa
Specificity	The mouse monoclonal antibody BA5b recognizes an extracellular epitope of CD26, a 110 kDa type II transmembrane glycoprotein, which is a peptidase expressed on mature thymocytes, T cells (especially activated), B cells, NK cells and macrophages.
Other names	DPP4, TP103, ADCP-2, ADABP
Workshop	HLDA VI: WS Code N-L078
Immunogen	A human T cell clone
Entrez Gene ID	1803
Gene name	DPP4
NCBI Full Gene Name	dipeptidyl peptidase 4
UniProt ID	P27487

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.