



	Technical Data Sheet	
Product	Anti-Hu CD42a Purified	
Cat. Number/Size	11-732-C025	0.025 mg
	11-732-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD42a	
Clone	GR-P	
Format	Purified	
Reactivity	Dog, Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody GR-P (also known as GRP-P) recognizes an extracellular epitope of CD42a (glycoprotein 9), a 22 kDa transmembrane protein constitutively expressed on megakaryocytes and platelets.	
Other names	GPIX, GP9	
Workshop	HLDA IV: WS Code P 35	
Immunogen	Human acute lymphoblastic leukemia cells	
Entrez Gene ID	2815	
Gene name	GP9	
NCBI Full Gene Name	glycoprotein IX platelet	
UniProt ID	P14770	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.