

Technical Data Sheet

Product

Anti-TAp63 Purified

Cat. Number/Size

11-890-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	TAp63
Clone	TAp63-4.1
Format	Purified
Reactivity	Mouse, Rat, Human
Application	WB, IHC-P, ICC, IP
Application details	Western blotting: Recommended dilution: 1 µg/ml
Isotype	Mouse IgG2a
Specificity	The mouse monoclonal antibody TAp63-4.1 recognizes TAp63; target epitope LSDPMW. This antibody does not bind to deltaNp63 isoform of p63.
Other names	AIS, KET, LMS, NBP, RHS, EEC3, OFC8, SHFM4
Immunogen	Recombinant human p63
Entrez Gene ID	8626
Gene name	TP63
NCBI Full Gene Name	tumor protein p63
UniProt ID	Q9H3D4

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.