

Technical Data Sheet

Product	Anti-AGPS Purified	
Cat. Number/Size	11-739-C025	0.025 mg
	11-739-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	AGPS
Clone	AGPS-03
Format	Purified
Reactivity	Human
Application	FC, WB
Isotype	Mouse IgG2a
Specificity	The mouse monoclonal antibody AGPS-03 recognizes AGPS (alkylglycerone phosphate synthase), an intracellular peroxisomal enzyme important for lipid biosynthesis.
Other names	alkylglycerone phosphate synthase, ADAS, ADPS, RCDP3, ADAP-S, ADHAPS, ALDHPSY
Immunogen	recombinant human AGPS (amino acids 158-384)
Entrez Gene ID	8540
Gene name	AGPS
NCBI Full Gene Name	alkylglycerone phosphate synthase
UniProt ID	O00116

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References www.exbio.cz

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.