

Technical Data Sheet

Product **Anti-Hu CD230 / Prion Purified**
Cat. Number/Size **11-604-C100** **0.1 mg**
For Research Use Only.
Not for use in diagnostic or therapeutic procedures.

Antigen	CD230 / Prion
Clone	EM-21
Format	Purified
Reactivity	Human
Application	WB
Application details	Western blotting: Recommended dilution: 0.5 µg/ml, non-reducing conditions are essential.
Isotype	Mouse IgG2a
Specificity	The mouse monoclonal antibody EM-21 recognizes human prion protein (PrP). Diglycosylated form of PrP has ~ 40 kDa, monoglycosylated form ~ 30 kDa, and nonglycosylated form ~ 19-21 kDa.
Other names	Major prion protein, PrP, PrP27-30, PrP33-35C, ASCR, CD230, PRNP, CJD, KURU, PRIP, PRPC
Immunogen	Recombinant human prion protein
Entrez Gene ID	5621
Gene name	PRNP
NCBI Full Gene Name	prion protein
UniProt ID	P04156

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.