

Technical Data Sheet

Product **Anti-HLA-DR Purified Low Endotoxin**  
Cat. Number/Size **12-690-C100** **0.1 mg**  
**For Research Use Only.**  
**Not for use in diagnostic or therapeutic procedures.**

Antigen HLA-DR  
Clone L243  
Format Purified Low Endotoxin  
Reactivity Dog, Non-human primates, Human  
Application FC (QC tested), IP, WB, IHC(P), IHC(F), ICC, FUNC  
Application details Functional application: Blocking.  
Flow cytometry: Recommended dilution: 1-4 µg/ml  
Isotype Mouse IgG2a kappa  
Specificity The mouse monoclonal antibody L243 recognizes specifically an extracellular epitope on HLA-DR molecules, both peptide-loaded and empty.  
Immunogen Human B lymphocytes

Concentration 1 mg/ml  
Preparation Purified by protein-A affinity chromatography.  
Formulation Phosphate buffered saline (PBS), pH 7.4  
Storage and handling Store at 2-8°C. Do not freeze.

Images and References [www.exbio.cz](http://www.exbio.cz)

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.