

Technical Data Sheet

Product	Anti-GRB2 Purified	
Cat. Number/Size	11-554-C025	0.025 mg
	11-554-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	GRB2
Clone	PAb (554)
Format	Purified
Reactivity	Human
Application	WB
Application details	Western blotting: Recommended dilution: 1-2 µg/ml, reducing conditions.
Isotype	Rabbit polyclonal
Specificity	The polyclonal antibody recognizes Grb2, a highly conserved approximately 26 kDa intracellular adaptor protein. The immunogen has 100% identity between man, cow, rat and mouse.
Other names	ASH, NCKAP2
Immunogen	Synthetic peptide corresponding to C-terminus of human Grb2 (amino acids 198-217).
Entrez Gene ID	2885
Gene name	GRB2
NCBI Full Gene Name	growth factor receptor bound protein 2
UniProt ID	P62993

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.