



	Technical Data Sheet	
Product	Anti-GFP Purified	
Cat. Number/Size	11-476-C025	0.025 mg
	11-476-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	GFP	
Clone	PAb (476)	
Format	Purified	
Reactivity	Tagged fusion proteins in all species	
Application	IP, WB, ICC	
Application details	Western blotting: Recommended dilution 0.5-1.5 µg/ml; positive control: transfected cells; negative control: non-transfected cells. Immunocytochemistry: Recommended dilution 1-3 µg/ml.	
Isotype	Rabbit polyclonal	
Specificity	The polyclonal antibody recognizes GFP, EGFP, EYFP fusion proteins in all species.	
Other names	Green fluorescent protein, EGFP, EYFP	
Immunogen	EGFP, a native full-length protein	
Gene name	GFP	
NCBI Full Gene Name	green fluorescent protein	
UniProt ID	P42212	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.