

Technical Data Sheet

Product	<b>Anti-SOCS3 Purified</b>	
Cat. Number/Size	<b>11-406-C025</b>	<b>0.025 mg</b>
	<b>11-406-C100</b>	<b>0.1 mg</b>

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

---

Antigen	SOCS3
Clone	SO1
Format	Purified
Reactivity	Human
Application	WB
Application details	Western blotting: Positive control: HeLa human cervix carcinoma cell line, reducing conditions.
Isotype	Mouse IgG2b
Specificity	The antibody SO1 reacts with SOCS3, an intracellular cytokine signaling inhibitor.
Other names	SOCS-3, CIS3, SSI3
Immunogen	Full-length SOCS3 protein.
Entrez Gene ID	9021
Gene name	SOCS3
NCBI Full Gene Name	suppressor of cytokine signaling 3
UniProt ID	Q6FI39

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.