

Technical Data Sheet

Product **Anti-Hu CD218a Purified**
Cat. Number/Size **11-887-C100** **0.1 mg**
For Research Use Only.
Not for use in diagnostic or therapeutic procedures.

Antigen	CD218a
Clone	H44
Format	Purified
Reactivity	Human
Application	WB, IHC-P, FC (QC tested), IP, IHC-F
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody H44 recognizes an extracellular epitope of CD218a, the alpha subunit of IL-18 receptor.
Other names	IL18R1, IL-18 R1, IL18-RA
Immunogen	Human CD218a expressed on a transfected NK0 cell line
Entrez Gene ID	8809
Gene name	IL18R1
NCBI Full Gene Name	interleukin 18 receptor 1
UniProt ID	Q13478

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.