

Technical Data Sheet

Product **Anti-Tb10.4 (M. tuberculosis) Purified**
Cat. Number/Size **11-186-C100 0.1 mg**
For Research Use Only.
Not for use in diagnostic or therapeutic procedures.

| | |
|---------------------|--|
| Antigen | Tb10.4 (M. tuberculosis) |
| Clone | PAb (186) |
| Format | Purified |
| Reactivity | Mycobacterium tuberculosis |
| Application | WB, ELISA |
| Isotype | Rabbit polyclonal |
| Specificity | The rabbit polyclonal antibody to Tb10.4 reacts with Mycobacterium tuberculosis protein Tb10.4 (EsxH). |
| Other names | Rv0288, EsxH |
| Immunogen | Recombinant Tb10.4 protein produced in E. coli |
| Entrez Gene ID | 886603 |
| Gene name | esxH |
| NCBI Full Gene Name | ESAT-6-like protein EsxH |
| UniProt ID | P9WVK3 |

| | |
|----------------------|---|
| Concentration | 1 mg/ml |
| Preparation | Purified by ligand affinity chromatography. |
| Formulation | Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide |
| Storage and handling | Store at 2-8°C. Do not freeze. |

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.