



Technical Data Sheet

Product	Anti-HIV-1 capsid protein p24 Purified	
Cat. Number/Size	11-327-C025	0.025 mg
	11-327-C100	0.1 mg
For Research Use Only.		
Not for use in diagnostic or therapeutic procedures.		

Antigen	HIV-1 capsid protein p24
Clone	ND1
Format	Purified
Reactivity	Human Immunodeficiency Virus
Application	WB, ICC, ELISA
Isotype	Mouse IgG1
Specificity	The antibody ND1 reacts with HIV-1 (Human Immunodeficiency Virus Type 1) capsid protein p24 (gp24). Specificity of antibody ND1 was verified by Western blotting in reducing conditions with Sanofi Pasteur NEW LAV-BLOT I kit. Epitope was reported to be located within the amino acids 51-110.
Immunogen	HIV-1 p24 C-terminal peptide.
Concentration	1 mg/ml
Preparation	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.
Images and References	www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.