



	Technical Data Sheet	
Product	Anti-Hu CD23 Purified	
Cat. Number/Size	11-689-C025	0.025 mg
	11-689-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD23	
Clone	EBVCS-5	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP, IHC-F	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody EBVCS-5 recognizes an extracellular epitope located in the stalk region of human low affinity IgE receptor (CD23) between the 37 and 25 kDa cleavage sites.	
Other names	FCER2, Fc-epsilon-RII, BLAST-2, FCE2, CLEC4J, IGEBF	
Immunogen	EBV-transformed human cells	
Entrez Gene ID	2208	
Gene name	FCER2	
NCBI Full Gene Name	low affinity immunoglobulin epsilon Fc receptor	
UniProt ID	P06734	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.