

Technical Data Sheet

Product	<b>Anti-Ms CD5 Purified</b>	
Cat. Number/Size	<b>11-748-C025</b>	<b>0.025 mg</b>
	<b>11-748-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	CD5
Clone	53-7.3
Format	Purified
Reactivity	Mouse
Application	FC (QC tested), IP, IHC(P), IHC(F)
Isotype	Rat IgG2a kappa
Specificity	The rat monoclonal antibody 53-7.3 recognizes an extracellular epitope of CD5, a 67kDa single-chain transmembrane glycoprotein expressed on mature T lymphocytes, most of thymocytes and B-1 lymphocytes.
Other names	T1, LEU1
Immunogen	murine thymus or spleen cells
Entrez Gene ID	12507
Gene name	CD5
NCBI Full Gene Name	CD5 molecule
UniProt ID	P13379

---

Concentration	1 mg/ml
Preparation	Purified from cell culture supernatant by protein-G affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.