

Technical Data Sheet

Product	<b>Anti-CTBP1 Purified</b>
Cat. Number/Size	<b>11-511-C100</b> <b>0.1 mg</b>
	<b>For Research Use Only.</b>
	<b>Not for use in diagnostic or therapeutic procedures.</b>

---

Antigen	CTBP1
Clone	PAb (511)
Format	Purified
Reactivity	Mouse, Human
Application	IP, WB, ICC
Isotype	Rabbit polyclonal
Specificity	The polyclonal antibody reacts with CtBP1, a 49 kDa transcriptional corepressor (intracellular antigen), that plays a critical role in regulating gene expression.
Other names	BARS, HADDTS
Immunogen	Bacterially expressed full-length human CtBP1
Entrez Gene ID	1487
Gene name	CTBP1
NCBI Full Gene Name	C-terminal binding protein 1
UniProt ID	Q13363

---

Concentration	1 mg/ml
Preparation	Purified from rabbit serum by affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.