

Technical Data Sheet

Product

Anti-HLA-ABCE PE

Cat. Number/Size

1P-813-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	HLA-ABCE
Clone	TP25.99SF
Format	PE
Reactivity	Human
Application	FC (QC tested)
Excitation laser	blue (488 nm)
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody TP25.99SF recognizes HLA-ABC and HLA-E molecules, but not HLA-G. It can be used for discrimination between HLA-G and other HLA-class I antigens.
Immunogen	IFN-gamma-treated human melanoma cells Colo 38

Concentration	0.1 mg/ml
Preparation	The purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. The conjugate is purified by size-exclusion chromatography.
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store in the dark at 2-8°C. Do not freeze. Avoid prolonged exposure to light. Do not use after expiration date stamped on vial label.

Images and References www.exbio.cz

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.