

Technical Data Sheet

Product **Anti-Hu CD199 Purified**
Cat. Number/Size **11-902-C100** **0.1 mg**
For Research Use Only.
Not for use in diagnostic or therapeutic procedures.

Antigen	CD199
Clone	C9Mab-1
Format	Purified
Reactivity	Human
Application	FC (QC tested)
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody C9Mab-1 recognizes an extracellular epitope of CD199 (Chemokine Receptor 9), a 7-transmembrane chemokine receptor.
Other names	GPR-9-6, GPR28, CCR9
Immunogen	CD199 transfected CHO cells
Entrez Gene ID	10803
Gene name	CCR9
NCBI Full Gene Name	C-C motif chemokine receptor 9
UniProt ID	P51686

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.