

Technical Data Sheet

Product	Anti-Cytokeratin 18 Purified	
Cat. Number/Size	11-110-C025	0.025 mg
	11-110-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	Cytokeratin 18
Clone	DA-7
Format	Purified
Reactivity	Human
Negative species	Dog, Mouse, Rat, Pig, Hamster, Sheep, Cow
Application	IP, WB, IHC(P), ICC, ELISA
Isotype	Mouse IgG1
Specificity	The antibody DA-7 reacts with Cytokeratin 18 (45 kDa). Cytokeratins are a member of intermediate filaments subfamily represented in epithelial tissues.
Other names	K18, CK18, CK-18, CYK18, KRT18
Immunogen	Human breast carcinoma cell line PMC-42.

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on vial label.

Images and References www.exbio.cz

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.