

Technical Data Sheet

Product	Anti-Ms CD8a APC	
Cat. Number/Size	1A-579-C025	0.025 mg
	1A-579-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD8a
Clone	53-6.7
Format	APC
Reactivity	Mouse
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-2 µg/ml.
Excitation laser	red (633 nm)
Isotype	Rat IgG2a kappa
Specificity	The rat monoclonal antibody 53-6.7 recognizes an extracellular epitope of murine CD8a (32-34 kDa; alpha chain of the CD8 antigen).
Other names	Ly-2, Ly-B, Ly-35, Lyt-2, BB154331
Immunogen	Murine spleen cells
Entrez Gene ID	12525
Gene name	CD8a
NCBI Full Gene Name	CD8a molecule
UniProt ID	P01731

Concentration	0.5 mg/ml
Preparation	Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.