



	Technical Data Sheet	
Product	Anti-Hu CD112 Purified	
Cat. Number/Size	11-777-C025	0.025 mg
	11-777-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD112	
Clone	R2.525	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP, IHC-F	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody R2.525 recognizes an extracellular epitope on CD112, a type I transmembrane glycoprotein expressed by myelomonocytic and megakaryocytic cells, and by CD34+ hematopoietic progenitors.	
Other names	NECTIN2, HVEB, PRR2, PVRL2, PVRR2	
Workshop	HCDM listed	
Immunogen	NIH/3T3 cells transfected with human Nectin-2	
Entrez Gene ID	5819	
Gene name	NECTIN2	
NCBI Full Gene Name	nectin cell adhesion molecule 2	
UniProt ID	Q92692	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.