



Technical Data Sheet

Product	Anti-HSV-2 glycoprotein G Purified	
Cat. Number/Size	11-334-C025	0.025 mg
	11-334-C100	0.1 mg
For Research Use Only.		
Not for use in diagnostic or therapeutic procedures.		

Antigen	HSV-2 glycoprotein G
Clone	T303
Format	Purified
Reactivity	Herpes Simplex Virus
Application	WB, ICC, ELISA
Isotype	Mouse IgG2a
Specificity	The antibody T303 recognizes gG antigen of HSV2 (Herpes Simplex Virus type 2), a dsDNA virus, member of Herpesviridae family.
Other names	HSVGG, HSV-2 gG
Immunogen	Extract of HSV-2 infected VERO (green monkey kidney) cells.
Entrez Gene ID	1487356
Gene name	US4
NCBI Full Gene Name	envelope glycoprotein G
UniProt ID	A7U885

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References	www.exbio.cz
-----------------------	--

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.