

Technical Data Sheet

Product

Anti-Hu CD79a Purified

Cat. Number/Size

11-530-C025

0.025 mg

11-530-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	CD79a
Clone	HM57
Format	Purified
Reactivity	Chicken, Rabbit, Opossum, Guinea pig, Horse, Cow, Pig, Mouse, Rat, Human
Application	FC (QC tested), IHC(P), IHC(F)
Isotype	Mouse IgG1
Specificity	The antibody HM57 interacts with CD79a (Ig alpha), a 40-45 kDa subunit of B cell antigen-specific receptor (BCR) and its early developmental forms.
Other names	BCR alpha, Ig-alpha, MB-1, IGA
Workshop	HLDA V: WS Code BC cB018; HLDA VI: WS Code BP 193; HLDA VI: WS Code BP 89; HLDA VI: WS Code B B103; HLDA VI: WS Code B CD79.4
Immunogen	Synthetic peptide corresponding to amino acids 202-216 of human CD79a

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on vial label.

Images and References www.exbio.cz

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.