

Technical Data Sheet

Product	Anti-UBN1 Purified	
Cat. Number/Size	11-565-C025	0.025 mg
	11-565-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	UBN1
Clone	UBN1-02
Format	Purified
Reactivity	Human
Application	IP, WB
Application details	Western blotting: Recommended dilution: 1-5 µg/ml; positive control: HeLa cell line.
Isotype	Mouse IgG1 kappa
Specificity	The antibody UBN1-02 recognizes N-terminal part of ubinuclein 1 (UBN1), a widely expressed nuclear and adhesion complex protein. Western blotting analysis reveals UBN1 as a 165 kDa band.
Other names	Ubinuclein 1, VT4
Immunogen	Recombinant protein corresponding to amino acids 1-190 of ubinuclein 1.
Entrez Gene ID	29855
Gene name	UBN1
NCBI Full Gene Name	ubinuclein 1
UniProt ID	Q6PIZ9

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.