

Technical Data Sheet

Product	<b>Anti-Hu CD158z Purified</b>
Cat. Number/Size	<b>11-945-C100</b> <b>0.1 mg</b>
	<b>For Research Use Only.</b>
	<b>Not for use in diagnostic or therapeutic procedures.</b>

---

Antigen	CD158z
Clone	CH21
Format	Purified
Reactivity	Human
Application	FC (QC tested), IP
Application details	Flow cytometry: Recommended dilution: 3-12 µg/ml
Isotype	Mouse IgG2a
Specificity	The mouse monoclonal antibody CH21 recognizes an extracellular epitope of human CD158z (KIR3DL3), a transmembrane glycoprotein of killer cell inhibitory receptor family.
Other names	KIR3DL3
Immunogen	human CD158z transfectants
Entrez Gene ID	115653
Gene name	KIR3DL3
NCBI Full Gene Name	killer cell immunoglobulin like receptor, three Ig
UniProt ID	Q8N743

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.