

Technical Data Sheet

Product	Anti-Hu TCR gamma/delta Purified Low Endotoxin	
Cat. Number/Size	12-147-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	TCR gamma/delta
Clone	B1
Format	Purified Low Endotoxin
Reactivity	Human, Non-human primates
Application	FC (QC tested), IHC(P), IHC(F), FUNC
Application details	Functional application: Blocking. Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody B1 (also known as B1.1) recognizes an extracellular epitope of TCR gamma/delta, the subtype of T cell receptor expressed mainly in epithelial tissues and at the sites of infection.
Other names	TCRG/D

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.