



Technical Data Sheet

|                                                             |                                                                                                                                                                                                   |                 |
|-------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| Product                                                     | <b>Anti-Cytokeratins Alexa Fluor® 488</b>                                                                                                                                                         |                 |
| Cat. Number/Size                                            | <b>A4-108-C025</b>                                                                                                                                                                                | <b>0.025 mg</b> |
|                                                             | <b>A4-108-C100</b>                                                                                                                                                                                | <b>0.1 mg</b>   |
| <b>For Research Use Only.</b>                               |                                                                                                                                                                                                   |                 |
| <b>Not for use in diagnostic or therapeutic procedures.</b> |                                                                                                                                                                                                   |                 |
| <hr/>                                                       |                                                                                                                                                                                                   |                 |
| Antigen                                                     | Cytokeratins                                                                                                                                                                                      |                 |
| Clone                                                       | C-11                                                                                                                                                                                              |                 |
| Format                                                      | Alexa Fluor® 488                                                                                                                                                                                  |                 |
| Reactivity                                                  | Mammalian                                                                                                                                                                                         |                 |
| Application                                                 | FC (QC tested), WB, ICC                                                                                                                                                                           |                 |
| Application details                                         | Flow cytometry: Recommended dilution: 1 µg/ml. Intracellular staining.                                                                                                                            |                 |
| Excitation laser                                            | blue (488 nm)                                                                                                                                                                                     |                 |
| Isotype                                                     | Mouse IgG1                                                                                                                                                                                        |                 |
| Specificity                                                 | The antibody C-11 reacts with cytokeratin peptides 4, 5, 6, 8, 10, 13, 18. Cytokeratins are members of intermediate filaments subfamily intracellular proteins represented in epithelial tissues. |                 |
| Other names                                                 | cytokeratin, CYK, CK, KRT                                                                                                                                                                         |                 |
| Immunogen                                                   | Keratin-enriched preparation from human epidermoid carcinoma cell line A431.                                                                                                                      |                 |
| <hr/>                                                       |                                                                                                                                                                                                   |                 |
| Concentration                                               | 0.1 mg/ml                                                                                                                                                                                         |                 |
| Preparation                                                 | Purified antibody is conjugated with Alexa Fluor® 488 NHS ester under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.            |                 |
| Formulation                                                 | Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide                                                                                                                           |                 |
| Storage and handling                                        | Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.                                                                                                                          |                 |
| <hr/>                                                       |                                                                                                                                                                                                   |                 |
| Images and References                                       | www.exbio.cz                                                                                                                                                                                      |                 |

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.

Alexa Fluor®, Pacific Blue™ and Pacific Orange™ are registered trademarks of Life Technologies Corporation.

This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not (1) use this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; or (c) manufacturing or quality assurance or quality control, and/or (2) sell or transfer this product or its components for resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5781 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com.