



	Technical Data Sheet	
Product	Anti-Hu SIGLEC10 Purified	
Cat. Number/Size	11-932-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	SIGLEC10	
Clone	5G6	
Format	Purified	
Reactivity	Human	
Application	IHC, ELISA, IP	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody 5G6 recognizes an extracellular epitope of human SIGLEC10, a sialic acid-binding lectin expressed on subsets of human leucocytes.	
Other names	PRO940, SGL2	
Immunogen	SIGLEC10 extracellular domain fused with human IgG1 Fc fragment	
Entrez Gene ID	89790	
Gene name	SIGLEC10	
NCBI Full Gene Name	sialic acid binding Ig like lectin 10	
UniProt ID	Q96LC7	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.