



Technical Data Sheet

Product	Anti-Blood Group A1B Culture supernatant	
Cat. Number/Size	99-807-L001	1.0 ml
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	Blood Group A1B	
Clone	HE-24	
Format	Culture supernatant	
Reactivity	Human	
Application	AGG	
Isotype	Mouse IgM	
Specificity	The antibody HE-24 distinguishes blood group A ₁ B from A ₂ B. The specificity of the antibody HE-24 was verified on >1000 samples of blood.	
Immunogen	Glycoprotein fraction isolated from saliva of secretors with blood group A.	
Preparation	Hybridoma culture supernatant concentrated by membrane ultrafiltration.	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.