



	Technical Data Sheet	
Product	Anti-Hu CD266 PE	
Cat. Number/Size	1P-755-C025	0.025 mg
	1P-755-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD266	
Clone	ITEM-4	
Format	PE	
Reactivity	Mouse, Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 2-4 µg/ml.	
Excitation laser	blue (488 nm)	
Isotype	Mouse IgG2b	
Specificity	The mouse monoclonal antibody ITEM-4 recognizes an extracellular epitope of CD266 / TWEAK R, a TNFR superfamily receptor for CD255 / TWEAK, a TNF-like weak inducer of apoptosis.	
Other names	TWEAKR, FN14, TNFRSF12A	
Immunogen	human CD266-transfected P815 cells	
Entrez Gene ID	51330	
Gene name	TNFRSF12A	
NCBI Full Gene Name	TNF receptor superfamily member 12A	
UniProt ID	Q9NP84	
Concentration	0.1 mg/ml	
Preparation	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.