

Technical Data Sheet

Product

Anti-TPX2 Purified

Cat. Number/Size

11-818-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	TPX2
Clone	TPX2-01
Format	Purified
Reactivity	Rat, Mouse, Human
Application	IP, WB, IHC(P), ICC, ELISA
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody TPX2-01 recognizes the epitope EPFVPKKEKKS (aa 636-646 in human) of TPX2, a microtubule-associated critical regulator of mitosis, which is overexpressed in many types of tumors and is a marker of worse prognosis in various cancers.
Other names	DIL2, HCTP4, REPP86
Immunogen	Human TPX2

Concentration	1 mg/ml
Preparation	Purified from cell culture supernatant by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on vial label.

Images and References www.exbio.cz

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.