



	Technical Data Sheet	
Product	Anti-Hu CD263 Purified Azide Free	
Cat. Number/Size	10-518-C025	0.025 mg
	10-518-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD263	
Clone	TRAIL-R3-02	
Format	Purified Azide Free	
Reactivity	Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The antibody TRAIL-R3-02 reacts with TRAIL-R3, a 35 kDa GPI-anchored extracellular membrane protein expressed mainly on neutrophils.	
Other names	TNFRSF10C, DcR1, TRAIL-R3, TRID, LIT	
Workshop	HLDA IX	
Immunogen	TRAIL-R3-hlgGhc fusion protein	
Entrez Gene ID	8794	
Gene name	TNFRSF10C	
NCBI Full Gene Name	TNF receptor superfamily member 10c	
UniProt ID	O14798	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.