



| | | |
|-----------------------|---|---------------|
| | Technical Data Sheet | |
| Product | Anti-PAX5 Purified | |
| Cat. Number/Size | 11-952-C100 | 0.1 mg |
| | For Research Use Only. | |
| | Not for use in diagnostic or therapeutic procedures. | |
| Antigen | PAX5 | |
| Clone | 1H9 | |
| Format | Purified | |
| Reactivity | Mouse, Human | |
| Application | WB (reported), ICC (reported), IP (reported), FC-IC (QC tested), IHC-P (reported) | |
| Application details | Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining. | |
| Isotype | Rat IgG2a kappa | |
| Specificity | The rat monoclonal antibody 1H9 recognizes transcription factor PAX5 (intracellular antigen). | |
| Other names | ALL3, BSAP | |
| Immunogen | GST-tagged recombinant mouse PAX5 | |
| Entrez Gene ID | 5079 | |
| Gene name | PAX5 | |
| NCBI Full Gene Name | paired box 5 | |
| UniProt ID | Q02548 | |
| Concentration | 1 mg/ml | |
| Preparation | Purified by protein-A affinity chromatography. | |
| Formulation | Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide | |
| Storage and handling | Store at 2-8°C. Do not freeze. | |
| Images and References | www.exbio.cz | |

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.