



	Technical Data Sheet	
Product	Anti-Hu CD276 Purified	
Cat. Number/Size	11-963-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD276	
Clone	7-517	
Format	Purified	
Reactivity	Human	
Application	ELISA (reported), FC (QC tested), IHC-F (reported)	
Application details	Flow cytometry: Recommended dilution: 0.5-4 µg/ml.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody 7-517 recognizes an extracellular epitope of human CD276 (B7H3) transmembrane glycoprotein.	
Other names	B7H3, B7RP-2	
Workshop	HLDA8: WS code DC ARN-80650	
Immunogen	monocyte-derived dendritic cells	
Entrez Gene ID	80381	
Gene name	CD276	
NCBI Full Gene Name	CD276 molecule	
UniProt ID	Q5ZPR3	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.