

Technical Data Sheet

Product	Anti-Cytokeratin 19 Biotin	
Cat. Number/Size	1B-120-C025	0.025 mg
	1B-120-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	Cytokeratin 19
Clone	A53-B/A2
Format	Biotin
Reactivity	Human
Application	IHC(P), FC (QC tested), IP, WB, ICC, ELISA
Isotype	Mouse IgG2a
Specificity	The mouse monoclonal antibody A53-B/A2 recognizes the rod domain of cytokeratin 19 (40 kDa), an intracellular antigen constituting intermediate cytoskeleton filaments. Cytokeratin 19 is not expressed in hepatocytes; it is often co-expressed with cytokeratin 7.
Other names	K19, CK19, CK-19, CYK19, KRT19, K1CS
Immunogen	MCF-7 human breast adenocarcinoma cell line
Entrez Gene ID	3880
Gene name	KRT19
NCBI Full Gene Name	keratin 19
UniProt ID	P08727

Concentration	1 mg/ml
Preparation	The purified antibody is conjugated with biotin-LC-NHS under optimum conditions. The reagent is free of unconjugated biotin.
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.