



	Technical Data Sheet	
Product	Recombinant human Interleukin-2	
Cat. Number/Size	EXB0040	0.05 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Format	Purified Azide Free	
Application	FUNC (verified)	
Description	Recombinant protein IL-2 is expressed in S2 cell (Drosophila). DNA sequence encoding 133 AAs was fused with Strep-tag at the N-terminus. A calculated molecular mass of recombinant protein is 16.7 kDa.	
Contents	Recombinant protein human Interleukin-2 - 1 vial (0.05 mg), lyophilized.	
Preparation	Purified by affinity chromatography. Lyophilized from a 0.2 µm filtered solution in Storage buffer (100 mM Tris, 150 mM NaCl, 0.01% Tween 20, pH 8.0).	
Storage and handling	Store at -20°C to -80°C. Reconstitute in sterile deionized water. Use reconstituted product immediately or aliquot for further storage at -20°C to -80°C.	
Images and References	www.exbio.cz	

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.