



	Technical Data Sheet	
Product	Anti-Rt CD8b FITC	
Cat. Number/Size	1F-656-C025	0.025 mg
	1F-656-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD8b	
Clone	341	
Format	FITC	
Reactivity	Rat	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Excitation laser	blue (488 nm)	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody 341 (also known as 34.1) recognizes rat CD8b, the 32-34 kDa beta chain of the CD8 coreceptor (extracellular epitope), expressed on T cell subsets and some other cell types, such as macrophages.	
Other names	Ly-3, Ly-C, Lyt-3	
Immunogen	CD8 positive Wistar rat splenic T cell hybridomas	
Entrez Gene ID	926	
Gene name	CD8b	
NCBI Full Gene Name	CD8b molecule	
UniProt ID	P10966	
Concentration	0.5 mg/ml	
Preparation	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.