

Technical Data Sheet

Product	Anti-Hu CD119 Purified	
Cat. Number/Size	11-875-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD119
Clone	GIR-208
Format	Purified
Reactivity	Human
Application	WB, FC (QC tested), IP
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody GIR-208 recognizes an epitope within extracellular part of CD119, a transmembrane glycoprotein of cytokine receptor family.
Other names	IFNGR1, IFNgammaRa, IMD27
Workshop	HLDA VI: C-110
Immunogen	CD119 purified from human placenta
Entrez Gene ID	3459
Gene name	IFNGR1
NCBI Full Gene Name	interferon gamma receptor 1
UniProt ID	P15260

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.