



	Technical Data Sheet	
Product	Anti-DLL4 Purified Low Endotoxin	
Cat. Number/Size	12-717-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	DLL4	
Clone	MHD4-46	
Format	Purified Low Endotoxin	
Reactivity	Human	
Application	FC (QC tested), FUNC	
Application details	Functional application: Blocking. Flow cytometry: Recommended dilution: 1-4 µg/ml	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody MHD4-46 recognizes the extracellular domain of DLL4 (Delta-like ligand 4), a type I transmembrane protein which plays an important role in vascular development.	
Other names	DL4, AOS6	
Immunogen	recombinant soluble human DLL4	
Entrez Gene ID	54567	
Gene name	DLL4	
NCBI Full Gene Name	delta like canonical Notch ligand 4	
UniProt ID	Q9NR61	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.