



	Technical Data Sheet	
Product	Anti-ABRA1 Purified	
Cat. Number/Size	11-619-C025	0.025 mg
	11-619-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	ABRA1	
Clone	ABRA1-01	
Format	Purified	
Reactivity	Human	
Application	IP, WB, FC-IC	
Application details	Flow cytometry: Recommended dilution: 1-2 µg/ml. Intracellular staining. Western blotting: Reducing conditions.	
Isotype	Mouse IgG2a kappa	
Specificity	The mouse monoclonal antibody ABRA1-01 recognizes human ABRA1 (Abraxas, CCDC98), an adaptor protein involved in DNA repair, which migrates as a 45 kDa band on PAAGE under reducing conditions.	
Other names	ABRA1, Abraxas 1, CCDC98, FAM175A	
Immunogen	Recombinant protein corresponding to amino acids 1-313 of ABRA1 with N-terminal His6 tag	
Entrez Gene ID	84142	
Gene name	ABRAXAS1	
NCBI Full Gene Name	abraxas 1, BRCA1 A complex subunit	
UniProt ID	Q6UWZ7	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.