



	Technical Data Sheet	
Product	Anti-SCIMP Purified	
Cat. Number/Size	11-704-C025	0.025 mg
	11-704-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	SCIMP	
Clone	NVL-07	
Format	Purified	
Reactivity	Human	
Application	FC-IC, IP, ICC, WB (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining. Western blotting: Recommended dilution: 1-2 µg/ml.	
Isotype	Mouse IgG2a	
Specificity	The mouse monoclonal antibody NVL-07 recognizes intracellular part of human transmembrane adaptor SCIMP. This protein of 17 kDa predicted Mw migrates as a 22 kDa band on SDS PAGE.	
Other names	NVL	
Immunogen	Recombinant intracellular part of human SCIMP	
Entrez Gene ID	388325	
Gene name	SCIMP	
NCBI Full Gene Name	SLP adaptor and CSK interacting membrane protein	
UniProt ID	Q6UWF3	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.