

Technical Data Sheet

Product **Anti-Hu Galectin-9 Purified**  
Cat. Number/Size **11-928-C100** **0.1 mg**  
**For Research Use Only.**  
**Not for use in diagnostic or therapeutic procedures.**

---

Antigen	Galectin-9
Clone	9M1-3
Format	Purified
Reactivity	Non-human primates, Human
Application	FC-IC (QC tested), IHC (reported), FC (reported)
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody 9M1-3 recognizes an epitope within C terminus of human galectin-9, a glycan-binding protein expressed mainly on activated Th cells and FoxP3+ Treg cells.
Other names	HUAT, LGals9A, LEG9
Immunogen	recombinant M-type splicing variant of human galectin 9
Entrez Gene ID	3965
Gene name	LGALS9
NCBI Full Gene Name	galectin 9
UniProt ID	O00182

---

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.