



	Technical Data Sheet	
Product	Anti-Hu CD93 Purified	
Cat. Number/Size	11-794-C025	0.025 mg
	11-794-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD93	
Clone	VIMD2	
Format	Purified	
Reactivity	Human	
Application	WB, ICC, FC (QC tested), IHC, IP	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody VIMD2 recognizes an extracellular epitope on CD93, an approximately 110-120 kDa glycoprotein expressed mainly on myeloid cells and endothelial cells.	
Other names	C1QR1, C1qRP, ECSM3, MXRA4, C1qR(P), dJ737E23.1	
Workshop	HLDA I: WS Code M 3; HLDA V: WS Code M17, BP231; HLDA VI: WS Code MR8	
Immunogen	KG1 cell line	
Entrez Gene ID	22918	
Gene name	CD93	
NCBI Full Gene Name	CD93 molecule	
UniProt ID	Q9NPY3	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.